Link to different application instructions with NRX® Heat Tape and thermoplastics – www.mediroyal.se/en/brochures-for-nrx-strap-Instructions

© Copyright MediRoyal Nordic AB · Sweden · 2024

26205DR · www.vajer.se

mediroyal

NRX® Heat Tape Thermoplastic

MediRoyal Nordic AB
Box 7052 · SE-192 07 Sollentuna · Sweden
Tel +46 8 506 766 00 · www.mediroyal.se

English – instructions

NRX® Heat Tape Thermoplastic Kit

General information · The thermoplastic kit is made from 1,6 mm micro perforated thermoplastic that has been specially modified to work together with the NRX® Heat Tape products. The thermoplastic can be laminated on the NRX® Strap, ARX® Strap or SRX® Strap together with NRX® Heat Tape by using an iron.

Contains: 20 pcs of strips 15 x 200 mm
3 pcs of mini sheets 140 mm x 200 mm

Recommended temperature · Regular Iron set on medium heat, two dots

Working instruction · You need a silicone baking mat as surface protection on your working area. You also need a regular iron that can be set on medium heat (two dots) and regular baking paper.

Choose the quality of the NRX® Heat Tape that you wish to reinforce with thermoplastic. The regular NRX® Heat Tape has no hook compatibility while the NRX® Loop Heat Tape is hook compatible. Make sure that the NRX® Heat Tape covers the thermoplastic part with an overlap of at least 3 mm. Apply the thermoplastic to the tape surface and position it over the NRX® Strap, ARX® Strap or SRX® Strap material.

Cover the complete construction with a baking paper. Apply the iron over the materials and wait 10-15 seconds depending on the size of the construction. Check the result and if needed, repeat the process. If you use the regular NRX® Heat Tape, the thermoplastic will bleed through, and the surface becomes a bit sticky. Always check that the NRX® Heat Tape has bonded well, if not, continue to heat with the iron and add a bit of pressure.

When the construction is soft, apply it on the patient and let it cool.

Application instructions · Different application instructions with NRX® Heat Tape and thermoplastics can be found on www.mediroyal.se/en/brochures-for-nrx-strap-Instructions

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