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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Functional Knee orthosis

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.Z2040					
Size	S	M	L	XL	XXL	
Thigh circ. taken 15 cm above mid-patella	39/47	47/53	53/60	60/67	67/75	
Length cm	40	40	40	40	40	
Colour	dark grey					

Specify right or left

- MAINTENANCE**
- ☒ Do not bleach
 - ☒ No chemical cleaning
 - ☒ Do not iron
 - ☒ Do not tumble-dry
 - ☒ Washing instructions:

• Frame: Hand wash, maximum temperature 30°C, with mild soap (the use of a wash cloth is advisable); air dry, do not expose to hot temperatures; better if dried with a cloth.
• Polsterung: wash the padding (attached wash buttons) and hand wash at 30°C or lower, with mild soap, rinse thoroughly. Air dry far from extreme temperatures.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- ACL, MCL, LCL, PCL instability and ruptures, combined instability
- Complex instabilities of the knee joint, traumatic or degenerative
- Pre/Post-surgery both rehabilitative or conservative
- Ligament protection during sports

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

Frame:

- Light-weight and resistant, made of a special aluminum alloy, adjustable (facilitated modeling on thigh and calf)
- New anatomical design for a perfect adherence
- Front socket for easy fit (no need to put the brace on from the foot)

Polystritisches KLICKLACK joint with an ultra-low profile, easy to adjust using wedges WITHOUT any tool:

- Extension: 0°, 10°, 20°, 30°, 40°
- Flexion: 45°, 60°, 75°, 90°
- Block: 0°, 10°, 20°, 30°

Fastening system:

- Locking buckles with safety button to wear the orthosis without using the Velcro® fastening system (adjusted only at the first application)
- Divided and reinforced Velcro® straps for better adaptation
- Soft velveteen straps in MTP padding, seamless and with no adhesives

Padding of the frame:

- Hypoallergenic, anti-migration, and multi-grip material
- Low profile clickclack to ensure safe positioning (removable for cleaning purposes)
- Extra-thickness along the calf for optimal comfort
- Pre-installed anatomic medial condyle pad to provide excellent stability
- Two sets of condyle pads of different thicknesses are included

PUTTING ON THE APPLIANCE

ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 Unfasten all Velcro® straps to their maximum length and open the four snap buckles.
- 2 Have the patient sit with the knee bent at approximately 30°. Place the brace, fully extended, on the leg, centering the joints on top of the patella (fig. A).
- NOTE: the joints will result positioned posterior to the center of the knee, but they will center when straps are closed.
- 3 Close snap buckle n. 1 and tighten the corresponding Velcro® strap, make sure the tibial grip fits the leg well (fig. B).
- 4 Repeat the operation in the same way for straps n.2 and n.3 (fig. C), making sure the thigh grip fits well on the leg.
- 5 Tighten strap n.4 to Velcro® and close buckle n.5 while tightening the corresponding strap to Velcro® (fig. D).
- 6 Make the patient stand up and tighten strap n.6 to Velcro® (fig. E).
- 7 Ask the patient to walk. Check that thigh and tibia grips, straps, and condylar supports are well positioned on the leg, without creating excessive pressure.
- 8 If necessary, adapt brace to the limb of the patient:
 - a. Shape the frame so as to perfectly follow the injured limb (fig. F).
 - b. The thicknesses in correspondence of the condyles can be modified: the brace has pre-mounted medial anatomical pads, and a low anatomical one laterally; there are two high and one low replacement pads, included in the package.

SPECIFIC POSTERIOR CRUCIATE LIGAMENT APPLICATIONS

- a. For the posterior cruciate ligament follow the sequence shown in fig. G. Do not follow previously mentioned sequence (points 2-6).
- Remain in a standing position to fasten strap n.2, once fastened, sit down to complete the fastening sequence.

FLX-EXTENSION ADJUSTMENT

- 9 Sliding the two half-covers sideways (fig.H1), extract the central locking element (fig.H2) and insert/substitute the wedges to the desired degree (fig.H3)
- 10 Re-insert the central locking element (fig. I1) sliding the two half-covers towards the middle until you hear the safety lock "clack" (fig. I2)

NOTE: Insert the extension wedge (the one with one notch only) first and then the flexion wedge (the one with two notches) when adjusting the block to 0°/10°/20°/30°. If you encounter any difficulty in the insertion of the second wedge, scrape off any paint in the knee brace or file the notches of the wedge.

HOW TO REMOVE THE BUCKLES (IF NECESSARY)

- 11 Remove padding from the frame.
- 12 From the inside, unscrew the screw which corresponds to the buckle that has to be removed.
- 13 Remove stud and buckle.

PRODUCT REMOVAL

- 14 Loosen straps n.4 and n.6, without removing them.
- 15 Open snap buckles n.1, 2, 4 and 5.
- 16 Remove the brace.

SUBSEQUENT APPLICATIONS

- 1 Loosen straps n.4 and n.6, WITHOUT REMOVING THEM, and open the snap buckles.
- 2 Sitdown with knee flexed at approximately 30° and place the fully extended brace on the leg, centering the joints on the upper part of the patella (fig. L).
- 3 Close the snap buckle n.1 (fig.M), being careful the tibia grip adheres well to the leg.
- 4 Repeat the operation starting with buckle n.2 followed by n.3 (fig.N), being careful the thigh grip adheres well to the leg.
- 5 Tighten strap n.4 to Velcro® and close buckle n.5 (fig.O).
- 6 Stand up and tighten strap n.6 to Velcro® (fig.P).

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