GB **English – instructions**
Plantar Fasciitis Brace PF8

General information - Read the instruction carefully before you start to use the product. Persons with diabetes, circulatory problems or a sensitive skin should consult their doctor before use. If you are allergic, check the material content before the use.

The PF8 Designs

PF8 Standard - The PF8 Standard is designed for normal every day use, either in combination with shoes or bare foot. It can be worn either daytime or nighttime.

PF8 Active - The PF8 Active is designed for persons with higher activity mostly in combination with shoes. Due to the higher compression it should only be used during activity.

Indications - Plantar fasciitis, Bone spur, Heel pain due to collapsed fat pad. **ⓘ**

Material - The Plantar Fasciitis Brace is made from a high tech compression material which provides graded compression to the foot and plantar area. The material content is 60% Polyamide, 20% Nylon and 20% Elastane.

Function - The different compression zones provides the foot with off-loading the plantar area, it support the transversal arch, reduces local swelling and are enhancing circulation and healing. It can be used both during the day as well as night-time.

Night time use - Some patients can get a pain relief by using the brace during the night. Make sure that you have used the brace during day time for at least a week before using it at night. If you experience loss of sensation, a tingling feeling in the toes, increased pain or swelling of the toes take it off and only continue to use it during day time.

Pre-application - To select the right size use the measurement chart for reference. Always try the size on for the best fit. The compression is high and the brace will feel tight at first.

Application - Put the brace on like a sock. Make sure that you have put it on right. The size indication should be under the foot. The area around the heel and under the foot has a higher compression that will provide support to the heel. The area that runs in front of the heel and around the foot helps to unload and lift the attachment of the plantar fascia to reduce pain.

Removing the brace - Slide it off like a sock, inside out. Then turn it. The brace should not be stored inside out as the compression zones might be damaged then.

Tips on Plantar Fasciitis - Plantar fasciitis is usually caused by overstrain of the plantar band that runs like a bow string on a bow. When the foot is loaded the plantar band expands. Inward rotation of the foot (pronation) causes and increased load and irritation of the plantar band. This might lead to a smaller inflammation in the start that eventually can develop into a much more painful and long lasting problem if not treated early with the correct measures. Many patients wrongly get the "bone spur" which is rarely the case. A bone spur is a calcification of the heel attachment of the plantar fascia and can be a result of plantar fasciitis. The best way to treat your plantar fasciitis is to know what causes the problem.

- **Check your biomechanical foot status.** If you pronate we recommend that you use a medial posting wedge or an insole like MR1500 MOW. That insole will support the foot medially and reduce irritation and load on the plantar band. **ⓘ**
- **Stretch.** Stretching of the plantar fascia and the tendons running just behind your medial malleolus can reduce pain. The muscles running on the inside of the ankle are starting as tendons in the foot and ending as central muscles in the calf. That's often why you might also feel a cramp or pain in your calf when you suffer from plantar fasciitis. A good stretch exercise is this one from our Brace Guide. Stretch is excellent to do in the morning or before any major exercise. Download the Brace Guide on www.mediroyal.se/en and follow the two stretch exercises for calf on page 56. **ⓘ**
- **Ice for pain reduction.** Pain reduction can be reached by very simple means like ice. To apply a smaller cold pack or ice to the arch or the heel can reduce the pain effectively, especially if you have been walking all day long. Ice should be applied after activity when no further activity should be done afterwards and for maximum 20 minutes.

Caring instructions and precautions - We recommend that the brace should be washed once a week to avoid skin irritation. Machine wash 40 degrees C are recommended then let it air dry. Patients with sensitive skin, diabetes or poor circulation should check with their doctor before using. Should your problem become worse after using the brace, stop using and consult a doctor or a physical therapist. Patients with allergies should check the material content before use.

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© Copyright MediRoyal Nordic AB, Sweden · 2016 · 24919DR · www.vajer.se