



IMMO+ NS

Product Options

Sheet Thickness	Sheet Size	Perforation
1,6 mm	60 x 43 cm	Micro, Mini, Multi, Opti
2,0 mm	60 x 44 cm	Micro, Mini, Multi, Opti
2,5 mm	60 x 45 cm / 60 x 90 cm	Mini, Multi, Opti
3,2 mm	60 x 45 cm / 60 x 90 cm	Mini, Multi, Opti
4,0 mm	60 x 45 cm / 60 x 90 cm	Mini, Opti

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Low-Temperature thermoplastic material with non-sticky coating on both sides for the fabrication of custom-molded rigid splints, orthoses, and adaptive equipment.

Processing method

The best method for processing the low-temperature thermoplastic is to place the material into a temperature-controlled water bath. Leave the material in water heated to 70°C for +/- 1 minute. Depending on the thickness and perforation type of the material, heating time can be extended. If several sheets are placed in a water bath at the same time, it is recommended to place wax paper between the sheets to prevent the material from adhering.

- **Precautions:** Application of the material onto the patient must be performed only by a trained healthcare professional.
- **Warning:** The device is intended for single patient use only.

Instructions for use

After the material has been softened in the warm water as detailed above, lift the sheet out of the water bath and place it immediately on a flat, clean surface. Excess moisture should be removed by gently blotting the surface of the material with dry towel or cotton cloth. Before handling the material, hands should be dipped into cold water to prevent imprints, for ease of molding, and to facilitate hardening of the material.

After removing from the water bath, the material will cool off and the optimal temperature for application and molding is about 40°C. For complex applications that require longer periods of time, material that is warmer than 40°C may be needed. In such cases, appropriate skin protection (e.g., application of padding stockinet) is required. **Warning: Do not** apply the material directly onto open wounds, macerated or sensitive skin. Appropriate skin protection must be applied.

- **Cutting:** The material can be cut when hot or cold using scissors or a box cutter knife.
- **Bonding:** Non sticky coating on both sides prevents accidental self-bonding after heating. To bond the material, remove protective coating at desired locations by scraping the surface gently with a sharp object. By pressing firmly together the uncoated hot surfaces of the material, a permanent bond may be made.
- **Fastening:** The use of self-adhesive hook/loop fastening system is recommended.
- **Edge finishing:** Rough edges may be polished or levelled by hand after placing the edge of the finished product into hot water (70°C) for 5-10 seconds to soften.

For reheating and remolding the use of a heat gun can be recommended. **Precaution:** Do not spot-heat the material over 180°C.

Care and cleaning

Avoid prolonged exposure to ultra violet light by storing the material in a cardboard box. Splints made from the material may be cleaned, when necessary, with soap and lukewarm water. Dry thoroughly before reapplication.

Notice: Any serious incident that has occurred in relation to the medical device should be reported to the manufacturer and the competent authority of the Member State in which the user or patient is established.