mediroyal

Tonus Splint

TONUS 

mediroyal

CE

Size Chart

Pediatric 45–90 mm	Standard 90–130 mm	Extra long 90–130 mm	
-----------------------	-----------------------	-------------------------	---



Take the width measurement over the MCP-joints

English	10	GB
Svenska	12	SE
Suomi	15	FI
Deutsch	18	DE
Français	21	FR
Italiano	24	IT
Español	27	ES
Português	30	PT
Nederlands	33	NL
Eesti	36	ET
Lietuvos	39	LT
Latvijas	42	LV
Slovenski	45	SL
Polski	47	PL
Русский	50	RU
Chinese Traditional	54	CH
Chinese Mandarin	58	CM
Japanese	61	JP
Korean	66	KO
Malay	71	MY
Arabic	74	AR
Persian	78	PE


NRX[®]
STRAP
mediroyal

TON812 · Finger Separating Strap for TONUS®

Instruction

Mediroyal Tonus Splint

Patients with neurological disorders often suffer from spasticity. The involuntary muscle tension often changes rapidly and can be affected by external factors, both physiological as well as psychological like anxiety, mood, lifestyle and pain. Poor hand positioning and instability of the MCP-arch can also trigger more spasticity. The treatment philosophy of the Mediroyal Tonus splint is to provide a good positioning of the hand in combination with providing an individual dynamic wrist extension with the inflatable bladder. The inflation provides a slight pressure under the MCP-joints to support the arch and to extend the joints which helps to extend both the PIP- and DIP joints. The aim is to provide constant and even stretch over the flexor tendons in order to make the stretch reflex less hypersensitive. Other functions are to maintain passive wrist extension, support weak muscles, stretch contractures and pressure prevention. The upper extremity is a chain of joints starting from the shoulder all the way down to the DIP-joint of the fingers. It's important to assess all joints of the patient in order to get the maximum effect of the treatment. Combinations of different orthotic solutions are often the most efficient treatment. Under point 11 you can find some more information on how to combine other types of orthoses in the treatment.

1. General information · Read the instruction carefully. The product should be fitted by a medical professional like an occupational therapist, physical therapist, orthopedic technician or a medical doctor. Patients with diabetes, circulatory problems or sensitive skin should consult their doctor before use. If the patient is allergic, check the material content on the textile label.

2. Indications · The product has been designed to provide tonus reduction and positioning of the hand for patients with spasticity like CVA/Stroke, Multiple sclerosis or traumatic spinal cord injuries.

3. Contraindications · Patients that might increase spasticity by volar stimulation might not be able to use the Tonus splint.

4. Pre Application · To select the right size, take the width over the MCP-joints. If the patient has a very high tonus, it might be beneficial to try to inhibit the spasticity before application. Correct positioning of the shoulder, external rotation of the upper limb and a gentle extension and stretch of the wrist and fingers might make application easier.

Evaluate the maximal extension that can be reached in the wrist in combination with the positioning of the fingers. Maximal extension usually triggers more spasticity that's why it's important to find the sub maximal position for each individual patient. The more extension that is reached in the wrist, the less extension can be reached in the fingers. This is caused by the shortening of the flexor tendons. The whole length of the flexor tendons have to be distributed over the complete chain of joints from the lower arm, wrist to the DIP-joints. It's therefore a good idea to slightly flex the wrist in order to be able to extend the fingers. To prepare for application, start with deflating the bladder completely with the black tip on the pump. Then remove the wool cover. Bend the platform in the wrist area, just behind the adapter of the bladder to allow flexion. In order to prevent pressure from the end on the splint against the lower arm, bend the end slightly.

5. Ulnar and radial deviation · If your Tonus-splint is fitted with a rotation adapter under the bladder you can adjust ulnar and radial deviation. The rotation adapter is easily identified by the golden screws on the front part of the platform. Use the screwdriver included in the package. Unscrew the golden screws slightly in order to rotate the adapter with the bladder. Don't adjust the black center screw. Find the most suitable ulnar/radial deviation and then tighten the golden screws in order to fixate the adapter. If further adjustment is needed, simply open the golden screws and change the angle.

6. Application · Before inserting the platform in the wool cover again, make sure that the rubber bumper around the adapter is in position. Then slide the bladder and platform into the cover.

Make sure to apply the splint carefully on the patient in order not to trigger more spasticity. Always pay attention to the position of the thumb so it won't be hyper extended. The easiest way is to hold the splint against the arm, controlling the position of the thumb and then allowing the fingers to grip around the deflated bladder. Secure the first strap that runs around the wrist. Make sure that the hook part on the inside of the strap attaches to the platform surface. Then close the strap around the wrist. Do not over tighten. Then close the second strap around the lower arm and check that the proximal end doesn't cut into the lower arm.

7. Inflating the bladder · The cover has a yellow marked opening on the inside of the palmar part. Through this hole you can access the valve. Make sure to use the clear tip of the pump when inflating. Insert the tip into the valve and make a quarter of a turn in order to open the safety seal in the valve. Pump once to inflate the bladder. Pay attention the MCP joints and the arch. Pump another 1-2 times more in order to inflate the bladder further. Once you have pumped 2-3 times the bladder has reached a more circular shape, allowing support to the MCP arch and a dynamic extension to the wrist. The dynamic extension means that the splint platform no longer will have a tight fit over the volar part of the wrist. There is usually no need for further inflation. Wait for 3-5 minutes and then carefully try the extension of each finger. The full effect of the splint is usually reached within 10-15 minutes. At that time you should be able to extend at least one or more of Dig II-V. Now secure the strap around the MCP joints.

The effect might vary depending on the individual patient. If you haven't reached the full effect, try to inflate the bladder further.

8. Time of usage · The best effect is reached when the splint is worn for 1-2 hours each time. It's better to divide the time of usage to 1-2 hours, two to three times a day rather than using the splint for 4-6 hours consecutively. The full effect of the splint can only be reached after 3-6 months of constant use. Try to fit the splint sessions into the daily routines of the patient, usually in between the meals and daily activities. The Tonus-splint is excellent to use when the patient performs activities that usually trigger spasticity like gait training.

9. Effect · The effect of the Tonus-splint is to provide constant and even tension over the flexor tendons in order to make the spasticity reflex less hyper sensitive. This will gradually lead to that the spasticity is triggered less easily but will require a continuous treatment with the splint for 3-6 months.

The Tonus-splint will also provide a better hand positioning to maintain passive range of motion along with pressure prevention.

10. Tonus splint and the use with Botulinum toxin injections · Patients that are undergoing treatment with upper extremity botulinum injections can successfully use the Tonus-splint. The Tonus-splint can both be used prior to the injection to decrease spasticity as well during the complete treatment with the injections. After the injection has been done the Tonus splint will help to position the hand, especially after manual stretch exercise, in order to maintain the passive range of motion that has been gained during the therapy. The splint can be adjusted gradually in extension for the improvement that has been gained. The recommendation is to use the Tonus splint during the injections at daytime and to use a more stable resting hand splint in thermoplastic during the night. The Tonus splint can also be used to bridge the gap between injections in order not to lose the increased passive range of motion that has been gained during the injection therapy.

11. Combination with Erix Three shoulder and NRX® Strap · In order to provide more individual extension of the fingers, the NRX® Strap can be used. You can treat one or two fingers a time. For one finger application, use the 30 mm width and round the edge on one side. Cut a hole in the material, about 5 mm from the edge. Apply the strap over the finger and make sure

GB**SE**

that the position will be close to the MCP-joint. Make sure that the hole is not too tight as it might cause edema. Adjust the length of the strap in order to get the extension effect needed and secure the strap with a hook part onto the wrist strap of the Tonus-splint. If you want to combine two fingers at the same time, use the 50 mm width and cut two holes. Apply like the one finger application above.

The shoulder position is very important to reduce spasticity. Many patients usually sublux their shoulder anteriorly. ErixThree Neuro shoulder can assist in creating a posterior safe position of the shoulder joint. The ErixThree shoulder can also be combined with the NRX® Strap to create a supination strap that can be connected to the Tonus-splint.

12. Finger separation · If there is a need to separate the fingers, there is a special finger separation strap that has been developed for Tonus. The TON812 has four individual loops attached into a foam strap. Apply the loops onto each finger and let the strap run from the dorsal side to the volar part of the splint. Attach the hook parts to the platform or the foam straps on the splint.

13. Removal · To remove the splint, open the MCP strap on the splint. Use the black end of the pump and insert that into the valve. Deflate the bladder completely and remove the pump. Open the straps around the wrist and lower arm and remove the splint. Check the skin for pressure marks after the use.

14. Material content · If the patient is allergic, check the material content on the labels of the platform and cover before use.

15. Caring instructions for the cover · The cover is made from a special composition of wool fibers that enhance ventilation and moisture transportation and also provides pressure absorption. The cover can be machine washed in 60 degrees C using the enclosed laundry bag. The cover should not be tumbled, it should air dry. Wool is a natural fiber; therefore the cover might release lint or fibers during the first washings. Once the cover has dried after washing, you can use a hair brush to make the surface soft again. We recommend that the cover should be washed at least once a week. Replacement covers can be ordered separately.

16. Caring instructions and precautions · The platform can be wiped with a mild soap solution or surface alcohol for disinfection. We do not recommend washing the platform along with the cover in a washing machine.

Patients with sensitive skin, diabetes or poor circulation should be observed frequently for pressure marks. Patients with allergies should check the material content on the label before use.

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Mediroyal – Functional Swedish Design

Product development is a combination of ideas, material, design and know-how. At Mediroyal we have more than 20 years of experience with developing products that provides both the patient and the therapist with enhanced functionality, very day. Swedish design is unique in many ways. Quality, performance, functionality and design are all features that make a difference. That's why we are proud to label our products with "Functional Swedish Design". Because orthoses are our passion.

Mediroyal Nordic AB · Box 7052 · Staffans väg 6 B · SE-192 07 Sollentuna · Sweden
Tel +46 8 506 766 00 · info@mediroyal.se · www.mediroyal.se