General information · Read the instruction carefully. Patients with diabetes, circulatory problems or sensitive skin should consult their doctor before use. If the patient is allergic, check the material content on the white textile label.

Indications · The product has been designed to be used during the night time to stretch the fascia under the foot. That treatment can provide reduction in pain during the day when suffering from plantar fascitis and also bone spur.

Contraindications · Patients with circulatory problems should consult their doctor before use as the sock might cause swelling during the night.

Pre Application · To select the right size use the measurement chart. Make sure that the sock is not too tight as it might cause swelling.

Application · Open the hook strap on the upper part of the sock as well as on the toe tip. Pull the sock up over the lower leg. Position it just above the calf muscle, below the knee fold. Close the upper strap around the leg, don't over tighten as it might cause circulatory problems. Pull on the strap attached on to the toe tip. Apply the strap through the ring on the upper front side and attach. Adjust the tension of the strap so the ankle reaches at least 90 degrees, preferable 80 degrees but the most important is that the toes will be extended up towards the knee joint. That's when you reach the most effective stretch of the plantar fascia. It's important to understand that the stretch effect should be considered over a longer period of time. A night splinting of the plantar fascia can a long time, sometimes up to three months.

Removal · Open the strap attached to the ring in the front. Then open the upper strap around the calf. Pull off the sock and make sure to close all the straps when the sock is not in use. That will prevent the hook from hooking onto the surface of the sock.

Caring instructions and precautions · We recommend that the brace should be washed once a week. Follow the washing instructions on the white label, close the straps and always use a laundry bag.

Be careful if you will get up and walk during the night with the sock on as it can be slippery against some surfaces. There is an anti-slip surface on the tip of the sock.

Patients with sensitive skin, diabetes or poor circulation should be observed. Patients with allergies should check the material content on the label on the product.

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© Copyright MediRoyal Nordic AB, Sweden - 2016 - 25056DR - www.vajer.se