

Orthoses are our passion

NRX® STRAP KITS

The NRX® Strap kits are a new addition to the NRX® strap line of products. It's simply pre-cut applications that is ready to go within a minutes. They are made from the normal NRX® Strap with the same functionality.



NRX1004 · NRX® Strap Kit for Thumb Abduction

This kit can be used to support and stabilize the thumb in abduction. Can be used for neurological patients to prevent or treat light adduction contractures or be used as a functional CMC-1 support. The application enables the patient to get a functional position and it's even possible to walk with a crutch wearing it. Hook parts for secure application are included.

Indications: CMC-1 joint problems in the thumb, light adduction contractures.

Universal size



NRX1002 · NRX® Strap Kit for Supination


The pre-cut kit can be used to supinate the forearm on neurological patients. It's attached over the thumb and then the spiral goes up over the elbow joint in order to control and assist radius/ulna rotation. Hook parts for secure application is included.

Indications: Assisted supination or pronation on neurological patients or post-fracture treatments.

Universal size



Orthoses are our passion



NRX1005 · NRX® Strap Kit for Hip Rotation

This kit is designed to be used for patients with internal or external rotation of the hip. The pelvic strap and the thigh cuff provides attachment points for the rotation part that is then attach onto the cuff and strap. Patients with hip problems usually have a shortening of the hip adductor which requires several physical therapy sessions. The kit can be used both as a diagnostic tool to evaluate the effect on the patient as well as after the manual therapy treatment to maintain a constant dynamic stretch over the hip adductor. The external rotation strap will also provide individual pressure over the gluteus muscles which help to extend the hip joint. The hook parts on the pelvic belt and the thigh cuff are sewn on to the applications. One of the hook parts on the rotation strap is fixed in place and the other one is adjustable.

Indications: Short adductor muscles in the hip especially adductor longus and brevis, hip extension problems.

Universal size



NRX1001 · NRX® Strap Kit for Hammertoes

A pre-cut kit for hammertoe applications. Consists of a double loop and a single loop depending on the need along with an anchor strap that goes over the foot. The loop parts will help to flex the metatarsal joint and is then secured directly on the anchor strap that goes over the foot. Hook parts for secure application are included.

Indications: Hammertoe problems

Universal size



NRX1003 · NRX® Strap Kit for Hallux Valgus

This pre-cut kit is used to control abduction of the toe in Hallux Valgus patients. Realigning the joint will not only lead to better neurological muscle function, it can also assist in reactivation the muscular stability of the joint as well as offering off-loading of pain in the joint. The strap is attached over the big toe and individual tension is applied over the strap that is secured over the ankle. A transversal arch band is the applied to secure the application. A sock can be worn over the application, making it possible to wear actively in shoes with enough space. Hook parts for secure application are included.

Indications: Hallux Valgus

Universal size