Mediroyal Nordic AB · Box 7052 · Staffans väg 6 B · SE-192 07 Sollentuna · Sweden
Tel +46 8 506 766 00 · info@mediroyal.se · www.mediroyal.se

© Copyright MediRoyal Nordic AB, Sweden · 2019 · 25521DR · www.vajer.se

TFCC Elbow Unit

TFC105

Boa® is a registered trademark of BoaTechnology Inc.

English
Svenska

GB
SE

GB

English · Instructions

Mediroyal TFCC Elbow Extension Unit

This elbow orthosis has been designed to be connected to the TFCC-wrist orthosis and to provide more rotation control. It must be used together with one of the models TFC101, TFC102 or TFC103 and are not compatible with any other orthosis. The elbow unit is available in left or right design and in three sizes 1-2-3.

General information · Read the instruction carefully. The product should be fitted by a medical professional like an occupational therapist, physical therapist, orthopedic technician or a medical doctor. Patients with diabetes, circulatory problems or sensitive skin should consult their doctor before use. If the patient is allergic, check the material content in the instruction before use.

Indications · This elbow orthosis has been designed to be used together a TFCC wrist orthosis to provide more rotation stability after TFCC surgery. It only works in combination with the wrist orthosis and not separately.

Contraindications · The elbow extension unit should not be used for fractures in the elbow region. The orthosis should always be used together with one of the TFCC-wrist orthosis.

Assembly of the elbow unit · Note that the elbow orthosis should be applied with an overlap against the wrist orthosis. This might make the short TFC101 model too short. The overlap needs to be 2-3 cm to provide a functional stability.

- Start by mounting the Boa® system onto the elbow orthosis. The Boa® dial should be applied on top of the overlapping lid. The two hook parts that are applied onto the wire should be positioned on the lower part of the elbow unit. Separate the parts slightly as this provides more stability.
- The condyle wings on the elbow unit can be fitted with extra pads for more compression over the elbow joint. Position the pad and attach the hook part over the edge to keep it in place. For more pronation control it might be enough with a pad on the medial side only.

Application and adjustment · The elbow unit should overlap the wrist orthosis with at least 2-3 cm for optimal stability.

- Open the Boa® system by lifting the dial until you hear a click. Pull out the wire completely.
- Remove the Boa® system that is attached with hook on the overlapping lid. Open the sides of the elbow unit and either sliding in from behind the elbow joint or over the wrist. Position the elbow unit so the condyle wings cover the side of the elbow joint and that the overlapping part over the wrist orthosis is at least 2-3 cm.
- Attach the Boa® system onto the overlapping lid again. Make sure that the wire runs over the surface of the unit and not over bare skin.
- Press the Boa® dial and adjust the wire with a light compression. Once again check the position of the elbow unit and then adjust more compression until it feels stable. Readjust if needed.
- The compression over the elbow joint can be reinforced by using the enclosed 5 cm NRX® Strap and the hook parts. The strap can be applied to provide more stability over the elbow joint. Attach the first hook over the medial condyle. Pull the strap diagonal and attach on the lateral side or attach it straight over to the lateral condyle.

A strap can also be attached on the backside of the elbow unit to provide more compression. This can be an alternative application if the patient needs unlimited extension and flexion in the elbow joint.

Additional rotation control with NRX® Strap · To provide additional and more specific rotation control the NRX® Strap or the NRX® PLUS Strap can be used by applying it like a spiral over the elbow unit and wrist orthosis.

Time usage · After the surgery the orthosis should be worn both night- and daytime or according to the specific instructions from your medical practitioner. It's also your medical practitioner that decides for how long you will wear the orthosis. The orthosis should be removed on a daily basis for ventilation of the skin and hygiene according to the instructions provided by your medical practitioner. This is important to prevent skin- or pressure problems.

Material content · The elbow extension unit is made from thermoplastic polymer and natural cork, polyamide and polyester.

Caring instructions and precautions · The elbow unit can be wiped with a mild soap solution and let air dry. It shouldn't be machine washed.

Patients with sensitive skin, diabetes or poor circulation should be observed for pressure marks or swelling. Patients with allergies should observe the material content. Should discomfort, skin irritation, increased pain or reduced sensation occur during your treatment, immediately contact your medical practitioner.

For more information visit www.mediroyal.se

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