



mediroyal

Orthoses are our passion

Manex Lace Up Wrist & Thumb

– Stable, breathable orthosis for best comfort, with a tailored fit!

MR2272 · Manex Lace-Up Wrist & Thumb

The Manex Lace-Up Wrist & Thumb is a new model that combines both thumb- and wrist stability with a lace-up closure system for easy application. The model is made from a Ventus™ material with loop surface on the outside. This makes it possible to apply the NRX® hook on the lace-up system anywhere on the surface of the brace. The Ventus™ material has a micro fleece lining, is comfortable and breathable.

The lace-up system is designed with two NRX® hook closures that provide a low profile and a secure attachment. The hook closures can be attached anywhere on the surface and can provide different levels of support depending on where they are positioned. The brace has a fixed tongue that makes it possible to apply like a glove. The thumb grip is adjustable.

Indications: General thumb and wrist instability, De Quervain, inflammatory diseases or tendinitis.

Manex Lace Up Wrist & Thumb Splints

Measurement chart

MR2272 Manex Lace-Up Short Wrist & Thumb
– Available in left or right model

Size	XS	S	M	L	XL	XXL	
cm	13–14	14–16	16–18	18–20	20–22	21–23	





Orthoses are our passion



MR2294 CTA Splint

Manex Lace Up Wrist & Thumb Splints

The splint for the wrist is the new CTA splint with a lower angle. The lower part of the splint has a split ending which decrease the pressure when the splint is loaded in volar flexion. The splint can be adjusted by using hot water.



MR2292 Splint

The thumb splint is pre-shaped made from anodized aluminum that can be easily adjusted.

To provide more stability the MR2490 Extra Strap 50 mm can be added.



MR2490 Extra Strap