



Orthoses are our passion



Omox Shoulder Immo Universal

MR4442 · Omox Shoulder Immo Universal

The Omox Shoulder Immo is a shoulder immobilizer that can be used on either left or right shoulder and is universal in size. It is made of ventilated foam and consists of an arm sling part that has two adjustable straps with hook, that should be crisscrossed in the back and attached on the sling for a secure positioning.

The straps have metal buckles that can be used for fine adjustments once the immobilizer is fitted.

To further stabilize the immobilizer a modular strap can be attached on the sling and wrapped around the body to prevent external rotation of the arm and to keep it fixed close to the body.

Indications: Post-operative immobilization of the elbow or shoulder or after trauma.

Measurement chart:

MR4442 · Omox Shoulder Immo Universal

Size

Universal size

