



Orthoses are our passion

Boa® Fit System for ErixThree® Neuro Shoulder

These are two modular systems that can be attached on to the ErixThree® Neuro shoulder brace to provide individual suspension and dynamic stability. The two systems have been designed with the Boa® Fit System for individual adjustability. The systems are attached directly on to the loop surface of the ErixThree® Neuro by the NRX® Hook that is applied on the tabs and the base of the systems. That allows you to individually position both the Boa® system as well as the tabs.

NRX945 2-point Boa® Fit System for ErixThree® Neuro

This system has two tabs that can be positioned individually. The Boa® system is integrated on a NRX® Hook reinforced with a carbon-nylon material. The excess hook can be trimmed if needed. The system has a return tab attached to create two individual wire loops where the NRX® tabs are attached to. The steel wire is coated with a polymer for protection.

The system can be attached over the trapezius for suspension or on the arm sleeve for dynamic stability. The tabs can be alternated to reach the optimal position for each individual patient.

NRX945 2-point Boa® Fit System for ErixThree® Neuro

Size

Universal size

NRX945

NRX946

See page 2

BOA®
Boa Fit System





BOA[®]
Boa Fit System

NRX946 3-point Boa[®] Fit System for ErixThree[®] Neuro

This system has three floating tabs that are positioned on the wire. The Boa[®] system is integrated on a NRX[®] Hook reinforcement with a carbon-nylon material. The excess hook can be trimmed if needed. The system is a single-wire system with floating tabs. The tabs can be positioned anywhere over the complete length of the wire. The steel wire is coated with a polymer for protection.

The system can be attached over the trapezius for suspension or on the arm sleeve for dynamic stability. The tabs can be positioned asymmetrically to reach the optimal position for each individual patient.

NRX946 3-point Boa[®] Fit System for ErixThree[®] Neuro

Size

Universal size



Orthoses are our passion

mediroyal

Vajer | wire | 2021 · 25801DR · ENG