

mediroyal

NRX[®]

BY MEDIROYAL

DOUBLE SOFT ELASTIC

The NRX[®] Double Soft Elastic is a hybrid strap between a normal elastic loop and an orthotic strap. It has a loop surface on both sides and can be cut without fraying. The NRX[®] Double Soft Elastic has a good stability yet conforms easily to the anatomy and the splint.

The strap has been developed to work together with the NRX[®] Low profile hook, both the self adhesive and the normal hook. This offers endless possibilities.

NRX[®] Double Soft Elastic is available in:

EU450080 NRX[®] Double Soft Elastic
80 mm width x 3,1 m

EU450110 NRX[®] Double Soft Elastic
110 mm width, 3,1 m

NRX[®] Low Profile Hook

EU445060 NRX[®] Low Profile Hook
Transparent Self-adhesive,
50 mm x 25 mm

EU445061 NRX[®] Low Profile Hook
Black with soft top surface,
50 mm width x 25 mm

