

MEDIROYAL PROTECTION MITTEN



HMT100

HMT100 MEDIROYAL PROTECTION MITTEN
Universally sized, wrist circumference 13-26 cm

The **Mediroyal protection mitten** has been designed to protect patients with disabilities from interfering with medical equipment like IV-tubes, tracheostomy tubes or catheters. It can also be used by patients who are prone to self-injury.

The Protection Mitten is made from a cotton material with an upper mesh surface. The cotton part is padded with cluster fiber that provides both comfort and cushioning for the palmar part of the hand. The palmar part has a slight dome shape which positions the hand in a reflex inhibitory position that can assist in mild spasticity. The upper part is made of a mesh material that provides good ventilation. The fingers have been separated by a seam that doesn't go all the way to the base of the finger to prevent pressure problems and skin irritation.

The closure of the Protection Mitten is made of a hook and loop strap that has a foam padded tongue to prevent pressure over the wrist. The hook and loop strap is secured through a plastic D-ring. In order to temporarily restrict the range of motion of the upper extremity there is a strap that can be fixated in the D-ring and then secured to the hospital bed. The Mediroyal Protection Mitten is universally sized and sold in pairs. It's machine washable in 40 degrees.