

MEDIROYAL OMOX ARM SLING



The Omox Arm Sling has been designed to provide an effective support and immobilization of the shoulder and arm after trauma.

It's made of a stable 3D material that both absorbs pressure and provides a good ventilation. The neck strap is adjustable and has a pad to absorb pressure. The sling has also a thumb strap to provide support to the thumb and wrist.

INDICATIONS

Soft tissue trauma to the shoulder, upper- or lower arm when immobilization is needed. Can be used in combination with a cast to off-load the additional weight and to provide immobilization.

MR1101-9 MEDIROYAL OMOX ARM SLING

The measurement is taken from the elbow tip to the MCP V joint.

Small	28-33 cm
Medium	33-40 cm
Large	40-51 cm