

# MEDIROYAL FOAM COLLAR



**The Mediroyal Foam Collar is a medium stable cervical collar made from medium density polyester foam.**

The collar has an anatomical shape for the chin to support the cheek. The back side is stable and has a closure with hook and loop. The collar is covered with two layers of cotton stockinette for maximum comfort.

#### Indications

To off-load the neck muscles or when light stabilization is needed.

**MR515 MEDIROYAL FOAM COLLAR GRAY  
MR516 MEDIROYAL FOAM COLLAR BLACK**

Height in the front	Gray	Black
6 cm	MR515-06	MR516-06
8 cm	MR515-08	MR516-08
10 cm	MR515-10	MR516-10
12 cm	MR515-12	MR516-12