

EDEMA GLOVE WITH SILICONE GRIP



TO IMPROVE THE GRIP FOR PATIENTS USING EDEMA GLOVES, Mediroyal has together with Hand Therapist Margareta Persson developed a glove with a grip improving pattern made of silicone. The pattern provides a high friction which leads to more functional grip for the patient. The silicone patterns are placed over each finger, thumb, thenar muscle and over vola. The glove has open fingers and is made in Lycra® with good compression and elastic

seams. The design provides excellent compression over the fingers and has no distal seam. The glove has a palmar pad included which can be used to improve compression over the MCP-joints or in the vola. Available in right and left version.

INDICATIONS Postoperative edema treatment or edema treatment after stroke.

906 EDEMA GLOVE WITH SILICONE GRIP, Right/Left

Size	Circumference around the MCP-joints
X-Small	178 mm
Small	203 mm
Medium	228 mm
Large	255 mm



DESIGN PATENTED BY MEDIROYAL NORDIC AB, SWEDEN