



mediroyal 

Orthoses are our passion

# Mediroyal Handy Wrist

– Short and versatile wrist support with Boa® Fit System



**BOA®**  
Boa Fit System

## MBA115 · Mediroyal Handy Wrist

The Mediroyal Handy Wrist is a short and versatile wrist support with a Boa® Fit System integrated to make application easy. The design is made with a ventilated fabric fitted with a profiled volar aluminium splint that can be adjusted individually. The unique design with a Lycra tongue makes application easy and the Boa® Fit System allows custom stability and compression.

The Mediroyal Handy Wrist is 15 cm long making it small, functional with a flexible stability and support. The Boa® Fit System makes application easy even for patients that lacks grip strength. The circular compression from the system provides a dynamic support that is easily adjusted.

**Indications:** Mild to moderate wrist instability, tendoniti

Size chart · Available in Left and Right version

### MBA115 · Mediroyal Handy Wrist

Size	S	M	L
Wrist circumference	12,5–16 cm	16–19 cm	19–23 cm



Boa® is a registered trademark of Boa Technology Inc.

Vajer |wire| 2022 - 26062DR - ENG