

THE ARMIS™ ANKLES BY MEDIROYAL



THE ARMIS™ LIGHT ANKLE BY MEDIROYAL



Flexible and low profile, yet strong and tough. That's a good description of our Armis™ Light Ankle.

It's made from a laminated air-mesh material that provides good ventilation and stability. The design has external reinforcements around the ankle to provide lateral and medial stability. The Boa® Closure System has been integrated on the back side to get a lower anterior profile. The system uses a polymer coated steel wire that provides stability and durability together with the rip-proof nylon guides in the front.

THE MOST VERSATILE

To provide more surface friction inside the brace we have added a semi-stretch sharkskin neoprene under the foot. The material has a loop function on the outside where the extra strap can be attached for individual support. The extra strap is designed to be attached over the Boa® system in the back to protect it from impact. The brace can also be used without the extra strap when a lighter support is needed. This combination makes the Armis™ Light Ankle one of the most versatile ankle braces on the market. Arm yourself with Armis™.



THE ARMIS™ ANKLE BY MEDIROYAL



High performance activities need functional and tough solutions.

That's why we have developed The Armis Ankle Brace by Mediroyal. It's made of high denier ballistic Nylon® with rip-stop function that can stand the toughest activities. The Armis™ Ankel is fitted with the Boa® Closure System for quick and easy custom support and stability.

Toughest on the market

The system uses a polymer coated steel wire that provides strength and durability together with rip-proof nylon and plastic wire guides. The design is seamless under the foot and has a flexible Achilles area, made from a breathable membrane. The lateral and medial stability comes from two plastic stays in the pockets of the brace. For individual stability the brace is fitted with two criss-cross straps that are anchored on the side of the brace. The combination with the Boa® Closure System, the plastic stays and the individual criss-cross strap makes the Armis™ Ankle one of the toughest on the market. Arm yourself with Armis™.



ARM YOURSELF WITH ARMIS™



MBA140 ARMIS™ LIGHT ANKLE

Size	Measurement taken around the malleolus
XX-Small	14-16 cm
X-Small	16-18 cm
Small	18-20 cm
Medium	20-22 cm
Large	22-24 cm
X-Large	24-26 cm

Boa® is a registered trademark of Boa Technology Inc.

MBA141 ARMIS™ ANKLE

Size	Measurement taken around the malleolus
X-Small	16-18 cm
Small	18-20 cm
Medium	20-22 cm
Large	22-24 cm
X-Large	24-28 cm
XX-Large	28-32 cm

Boa® is a registered trademark of Boa Technology Inc.

Indications Lateral or medial ankle instabilities, first to second degree acute ankle sprains or chronic ankle instability.

mediroyal 

www.mediroyal.se
info@mediroyal.se
facebook.com/Mediroyalnordic