



mediroyal

Orthoses are our passion



NRX

Thenax Opponent Thumb



MR2181 · NRX® Thenax Opponent Thumb


The NRX® Thenax Opponent is a custom moldable thumb for CMC-1 problems. It features an NRX® orthosis with a volar pocket for a thermoplastic part that can be custom molded to support the CMC-joint. The orthosis is made from NRX® material lined with Neptune™ textile for a high comfort. The pocket is made from 0,7 mm NRX® material where the thermoplastic is inserted. The thermoplastic part can be molded in the pocket by adding hot water or molded externally and then inserted into the pocket.

The orthosis has low profile hook for the thumb and wrist closure and an elastic hook strap that can be used to provide custom support over the CMC joint. The closure around the wrist features a D-ring that makes application easy.

Indications: CMC-1 instabilities, ligament injuries, inflammations, strains and sprains, rheumatoid arthritis, or degenerative diseases.

Measurement chart

MR2181 · NRX® Thenax Opponent Thumb · Available in left and right model

Size	XS	S	M	L	XL	
cm	13-14	14-16	16-18	18-20	20-22	

NRX® is a registered trademark by Mediroyal Nordic AB.

© Mediroyal Nordic AB 2021

Vajert/wirel 2021 · 25897DR · ENG