



mediroyal

Orthoses are our passion



## Origo Short PLUS

### MR2260P · Origo Short PLUS

Our popular Origo wrist is now available in a PLUS-size. The PLUS-size has been developed as a more tapered design with a wider proximal part to accommodate more conical arms. The design has no edging tape which makes the lower part of the construction more elastic. Compared to the regular Origo Basic model, the wrist measurement is the same and the proximal part will accommodate about 2-3 cm wider arms.

The article number for the PLUS size will be MR2260P. The design has four hook and loop closures, the upper one in 20 mm and the three others in 30 mm. The wrist splint is our dynamic plastic splint in the firm version, MR2297.

**Indications:** Wrist instabilities, post-operative after fractures, arthritis or tendonitis or for carpal tunnel syndrome.



Measurement chart · Available in left of right version

Size	S	M	L	XL
Wrist cm Ø	14–16	16–18	18–20	20–22
Proximal cm Ø	15–18,5	18–21,5	21–24,5	24–27,5

