

MEDIROYAL EPIARC™



EpiArc™ har desig och låg profil vid laterala epikondylitproblemen (tennisarm).

The EpiArc™ has been designed to provide an effective unloading and a low profile for lateral epicondylitis (tennis elbow). The large wing on the clasp provides a larger contact surface and has a silicone pad on the inside to provide an early signal from the extensor muscles over the lower arm. The silicone pad also provides a good grip to prevent the brace from migrating during activity. EpiArc can also be used for medial epicondylitis (golfers elbow). EpiArc™ has a lining made from NRX-material which provides a local heating effect. Heat provides more elasticity over muscle and can reduce pain further.

INDICATIONS

Lateral or medial epicondylitis.

MR2500 MEDIROYAL EPIARC™

Size	Measurement underarm
X-Small	21-24 cm
Small	24-25 cm
Medium	26-27 cm
Large	28-29 cm
X-Large	30-33 cm

