

# MEDIROYAL DORSUM FLEXI

MR4926



Dorsum Flexi is a new back support with a two-piece designed front closure. It's made from a ventilated elastic material that has been welded by using HF-technique.

This means that there are hardly any conventional seams in the construction. The benefits are a lower profile, higher comfort and a more durable construction. The two-piece front closure makes the height variable, which makes the model suitable for persons

with a wide hip and slim waist or with a larger abdomen. The lower closure has been reinforced to provide a better abdominal support, and has been fitted with handles to make the application easier. The upper straps have also been fitted with handles for easy application. The lumbar part has four splints, made of stainless steel, that have been laminated into a lumbar platform. The splints can easily be adjusted individually by bending but can't be taken out. There is a removable lumbar pad for additional support and comfort.

**INDICATIONS** Acute or chronic lumbago, lumbosacral instability, degenerative diseases, osteoporosis.

**MR4926 Mediroyal Dorsum Flexi**

Storlek	Circumference Waist
Small	75-83 cm
Medium	83-91 cm
Large	91-99 cm
X-Large	99-107 cm
XX-Large	107-115 cm
XXX-Large	115-125 cm

Height in the back: 27 cm.



# MEDIROYAL COTEX DORSUM STABLE

MR4928



A new flat knitted back-support with high compression and support, suitable for persons with high activity.

The front closure is reinforced to provide a good abdominal support. It's made from a HF welded material to provide a low profile and have handles to make application easy.

The reinforcement straps have been fitted with handles as well to make opening and closing

easier. The lumbar part has four splints, made of stainless steel, that have been laminated into a lumbar platform. The splints can easily be adjusted individually by bending but can't be taken out. The criss-cross strap is floating and can be applied with concentrated stability to the lumbosacral or lumbar area.

**INDICATIONS** Acute or chronic lumbago, lumbosacral instability, degenerative diseases, osteoporosis.

#### MR4928 Mediroyal CoTex Dorsum Stable

Size	Circumference Waist
Small	75-83 cm
Medium	83-91 cm
Large	91-99 cm
X-Large	99-107 cm
XX-Large	107-115 cm
XXX-Large	115-125 cm

Note that this model has a higher compression than normal back supports. This means that persons with less hand strength might have to use a larger size.

Height in the back: 26 cm  
Height in the front: 17 cm



mediroyal 

Mediroyal Nordic AB

Box 7052, 192 07 Sollentuna. Sweden

Visiting address: Staffans väg 6B. Sollentuna

Tel +46 8 506 766 00 Fax: +46 8 506 766 90

[www.mediroyal.se](http://www.mediroyal.se)

[info@mediroyal.se](mailto:info@mediroyal.se)

[facebook.com/MediroyalNordic](https://facebook.com/MediroyalNordic)