




mediroyal 

Orthoses are our passion

Extension 

EXTO[®]-15 Finger Ring - 15 degrees for Hyperextension

The new EXTO[®]-15 is a unique finger orthosis for hyperextension or hypermobility of finger joints. It can also be used for radial- or ulnar deviation to compensate and correct the deformity.



The EXTO[®]-15 has a 15 degree angled correction on each side of the ring in order to stop hyperextension earlier and more effectively. The length of the ring has been done a bit longer than the 2008 Multi Finger Splint, increasing the lever arm to provide more stability and support. The polymer composition will allow correction if necessary, by using a heat gun.

EXTO[®]-15 is available in sizes 1-14 with an asymmetrically graded length of each splint size. The measurement chart can be used as a guide to select the right size; however, an individual fitting is always recommended.

The EXTO[®]-15 is available as a single pack, a 3-pack and as an assortment pack with one size of each ring.

Indications: The product has been designed to support and protect the DIP and PIP finger joints against hyperextension, like swan neck deformity or hyper-mobility. It can also be used for radial- or ulnar deviation of the joints.



Size chart and available packages 

© Mediroyal Nordic AB 2024

Mediroyal Nordic AB · Box 7052 · 192 07 Sollentuna · Sweden · Tel +46 8 506 766 00
info@mediroyal.se · facebook.com/Mediroyalnordic · www.mediroyal.se

EXTO[®].15

EXTO[®]-15 Finger Ring - 15 degrees for Hyperextension

Available packages

Exto-15 Finger Ring Single pack

Size	Art.no.	Quantity
1	2015-01	One piece
2	2015-02	One piece
3	2015-03	One piece
4	2015-04	One piece
5	2015-05	One piece
6	2015-06	One piece
7	2015-07	One piece
8	2015-08	One piece
9	2015-09	One piece
10	2015-10	One piece
11	2015-11	One piece
12	2015-12	One piece
13	2015-13	One piece
14	2008-14	One piece


Exto-15 Finger Ring 3-pack

Size	Art.no.	Quantity
1	2015-01/03	3-pack
2	2015-02/03	3-pack
3	2015-03/03	3-pack
4	2015-04/03	3-pack
5	2015-05/03	3-pack
6	2015-06/03	3-pack
7	2015-07/03	3-pack
8	2015-08/03	3-pack
9	2015-09/03	3-pack
10	2015-10/03	3-pack
11	2015-11/03	3-pack
12	2015-12/03	3-pack
13	2015-13/03	3-pack
14	2015-14/03	3-pack

Exto-15 Finger Ring - Assortment Pack

Size	Art.no.	Quantity
1-14	2016	One piece of each size

One of the sides of the splint is a bit wider than the other part. This part is marked with a X-sign on the outside.



Size chart

Note that the measurement chart is a guide to select the right size. It does not take in consideration soft tissue swelling and edema. As most fingers have a conical shape the splint has one side that is slightly wider to compensate for this. However, we recommend an individual fitting for an optimal function.

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14
Circumference PIP joint mm	26	29	31	32	33	35	38	39	40	41	43	45	46	47
Length of splint mm	30	31,5	34	36	37	40,3	45	44	47,5	51	51,5	52	52,5	53

