




mediroyal 

Orthoses are our passion

Extension 

EXTO[®]-15 Finger Ring - 15 degrees for Hyperextension

The new EXTO[®]-15 is a unique finger orthosis for hyperextension or hypermobility of finger joints. It can also be used for radial- or ulnar deviation to compensate and correct the deformity.




The EXTO[®]-15 has a 15 degree angled correction on each side of the ring in order to stop hyperextension earlier and more effectively. The length of the ring has been done a bit longer than the 2008 Multi Finger Splint, increasing the lever arm to provide more stability and support. The polymer composition will allow correction if necessary, by using a heat gun.

EXTO[®]-15 is available in sizes 1-14 with an asymmetrically graded length of each splint size. The measurement chart can be used as a guide to select the right size; however, an individual fitting is always recommended.

The EXTO[®]-15 is available as a single pack, a 3-pack and as an assortment pack with one size of each ring.

Indications: The product has been designed to support and protect the DIP and PIP finger joints against hyperextension, like swan neck deformity or hyper-mobility. It can also be used for radial- or ulnar deviation of the joints.



Size chart and available packages 

© Mediroyal Nordic AB 2024

Mediroyal Nordic AB · Box 7052 · 192 07 Sollentuna · Sweden · Tel +46 8 506 766 00
info@mediroyal.se · facebook.com/Mediroyalnordic · www.mediroyal.se

EXTO[®].15

EXTO[®]-15 Finger Ring - 15 degrees for Hyperextension

Available packages

Exto-15 Finger Ring Single pack

Size	Art.no.	Quantity
1	2015-01	One piece
2	2015-02	One piece
3	2015-03	One piece
4	2015-04	One piece
5	2015-05	One piece
6	2015-06	One piece
7	2015-07	One piece
8	2015-08	One piece
9	2015-09	One piece
10	2015-10	One piece
11	2015-11	One piece
12	2015-12	One piece
13	2015-13	One piece
14	2008-14	One piece


Exto-15 Finger Ring 3-pack

Size	Art.no.	Quantity
1	2015-01/03	3-pack
2	2015-02/03	3-pack
3	2015-03/03	3-pack
4	2015-04/03	3-pack
5	2015-05/03	3-pack
6	2015-06/03	3-pack
7	2015-07/03	3-pack
8	2015-08/03	3-pack
9	2015-09/03	3-pack
10	2015-10/03	3-pack
11	2015-11/03	3-pack
12	2015-12/03	3-pack
13	2015-13/03	3-pack
14	2015-14/03	3-pack

Exto-15 Finger Ring - Assortment Pack

Size	Art.no.	Quantity
1-14	2016	One piece of each size

One of the sides of the splint is a bit wider than the other part. This part is marked with a X-sign on the outside.



Size chart

Note that the measurement chart is a guide to select the right size. It does not take in consideration soft tissue swelling and edema. As most fingers have a conical shape the splint has one side that is slightly wider to compensate for this. However, we recommend an individual fitting for an optimal function.

Size	1	2	3	4	5	6	7	8	9	10	11	12	13	14
Circumference PIP joint mm	42-43	44-46	46-50	50-55	55-58	58-60	60-62	62-65	65-68	68-71	71-74	74-77	77-81	81-84
Skenans längd mm	30	31,5	34	36	37	40,3	45	44	47,5	51	51,5	52	52,5	53

