



Orthoses are our passion

## Universal Wrist

UNI801 · UNI802

**– Bilateral wrist, universally sized**

The Universal wrist braces work for both left and right hand and are universally sized, covering the most common sizes.

The design is asymmetrical and combines the NRX-material with a foam laminate with high stability and a lycra sleeve to make application easy. Two splint pockets make it easy to switch the splint for right or left hand. The Universal wrist braces can be used as a base for NRX® Strap.

Universal wrist braces are available in two models.



UNI802



UNI801



### UNI801 Universal Wrist Short

Universal Wrist Short is 20 cm long and has three micro hook closures that can be attached directly to the brace. Supplied with the green CTA-splint MR2294.

**Indications:** CTA-Carpal tunnel syndrome, osteoarthritis, wrist instabilities, inflammations, sprains, rheumatoid arthritis, to off-load or hypermobility in the wrist.

### UNI802 Universal Wrist Long

Universal Wrist Long is 24 cm long and has two micro hook closures distal-proximal that can be attached directly to the brace. Around the wrist there is a 50 mm elastic loop strap that provides extra support. Supplied with the long stable splint MR2298.

**Indications:** After cast removal, wrist instability, inflammations, sprains or to off-load the wrist.



<input checked="" type="checkbox"/> UNI801 Universal Wrist Short	<input checked="" type="checkbox"/> UNI802 Universal Wrist Long	
Size	Universal size 12–24 cm	