

The Mediroyal logo is positioned at the top center of the image. It consists of the brand name 'mediroyal' in a white, lowercase, sans-serif font, followed by a stylized white graphic of three curved lines representing a wave or a spine. The logo is set against a dark blue rectangular background that also features a vertical yellow stripe on its right side.

mediroyal

Orthoses are our passion

Back Support Dorsum LOW

MR4924 • Back Support Dorsum LOW • 19 cm

The Dorsum LOW has been developed to provide extra support to the lumbar region in an extremely flexible and comfortable solution. The brace has a total height of 19 cm at the lumbar spine and only 14 cm at the front, giving the product a unique combination of flexibility and high comfort in one solution.

The back support is made of an elastic, thin, and breathable Lycra material with a thickness of 0.5 mm, ensuring optimal user comfort. Thanks to its unique construction, material, and low height, it can easily be worn underneath regular clothing.

Orthoses are our passion



In the lumbar area, four stays are positioned to locally increase support and relieve pressure. The stays are placed in open channels, two on the right and two on the left side of the spine. They can be removed, allowing the support to be individually adapted to the user's needs. The stays are made of flexible plastic that molds perfectly to the contours of the lower back.

In addition, the back support is equipped with an extra stabilizing X-strap across the lumbar region, making it easy and ergonomic to adjust the level of support as needed.

Recommended wash is 40° gentle cycle.

Close all hook-and-loop fasteners and remove the stays before washing.

After washing, lightly stretch and air-dry before use.

Measure waist circumference

MR4924 • Back Support Dorsum Low • 19 cm

Size	S	M	L	XL	XXL	XXXL
Cm	75-83	83-91	91-99	99-107	107-115	115-135