

CGO208

NRX[®] C-GO Thumb

A soft and elastic thumb orthosis in perforated CR-foam with a plush inside, that protects and off-loads the thumb. Designed without a splint to keep the thumb in a neutral resting position, making it particularly suitable for CMC- and MP-joint problems. The thumb support strap is made from Neptune[™] material and provides extra support for the CMC-joint, and the hook part under the strap secures the support function in the desired position.

The orthosis can be trimmed individually with a pair of scissors. The thumb part can be shortened or opened up if it is too tight. Do not cut over the seam.

CGO208 - Right / Left		Wrist circumference		
Size	1 (XS)	2 (S/M)	3 (L/XL)	4 (XXL)
cm	13-15	15-17	17-20	20-23

Indications: Osteoarthritis, CMC- or MCP-joint instability, inflammations, sprains, rheumatoid arthritis, off-loading or hyper mobility in the thumb or wrist.



CGO208

NRX[®] C-GO Kids

NRX[®] C-GO Kids is developed with the same proven design and functional principles as the NRX[®] C-GO Thumb, but is adapted to children's anatomical proportions.

There are two children's sizes.



CGO208 - Right/ Left		Wrist circumference
Size	00 (XXXS)	0 (XXS)
cm	9-11	11-13

