

DR127

Cervical Collar, Stable


A neck collar designed for long-term immobilization. Provides good stability in flexion, extension, lateral bending and axial rotation. The collar has flex zones in the front and the back to reduce pressure and increase comfort. The foam is latex-free and open cell to prevent skin breakdown. An extra set of pads is included.

Available in both pediatric and adult sizes. Available in five sizes.

Measure the circumference around the neck.

DR127					
Size	XS	S	M	L	XL
Circumference, cm	45-50	50-55	55-60	60-65	65-70
Height, cm	7	8	9,5	11,5	14

Indications: Cervical instabilities, fractures or whiplash injuries



ORTHOSES ARE OUR PASSION