



mediroyal 

Orthoses are our passion


**C-GO® Engineered design.
Functional support.**



NRX® C-GO Thumb

CGO208

A soft, elastic thumb orthosis made from perforated NRX® material with a plush inner lining for enhanced comfort. The outer surface features loop textile, allowing hook fasteners to attach securely at any point for a fully adjustable fit.



NRX®

CGO208 NRX® C-GO Thumb

The adjustable thumb sling strap is designed to protect and offload the CMC joint individually. Positioned over the thenar eminence, the sling includes a hook surface on the inside that attaches directly to the orthosis material, providing customized support, stabilization, and positioning of the CMC joint.

The orthosis features an easy-to-use two-piece hook closure system for simple application and secure fastening. For additional customization, the thumb section can be trimmed with scissors to achieve the desired length.

- Heat retention
- Adjustable support
- Pressure absorbing material

Indications: Osteoarthritis, CMC- or MCP-joint instability, inflammations, sprains, rheumatoid arthritis, after thumb arthroplasty or after surgery for thumb arthrosis, or hyper mobility in the thumb or wrist.

CGO208 – NRX® C-GO Thumb – Right / Left

Wrist circumference

Size	1 · XS	2 · S/M	3 · L/XL	4 · XXL	
cm	13–15	15–17	17–20	20–23	

Design protected in the EU under 15085791



Neptune™ C-GO Thumb

CG0108

Neptune™ is a low profile, breathable fabric engineered to provide a cooling sensation against the skin for enhanced all-day comfort. The Neptune™ C-GO Thumb features the same proven design as the NRX® C-GO, with a lower-profile construction and increased ventilation for improved comfort and wearability.

CG0108 Neptune™ C-GO Thumb

The adjustable thumb sling strap is specifically designed to support and offload the CMC joint. Positioned over the thenar eminence, the sling features an inner hook surface that attaches directly to the orthosis material, allowing for customized support, stabilization, and positioning of the CMC joint based on individual needs.

The orthosis also includes an easy-to-use two-piece hook closure system for quick application and secure fastening. For a more personalized fit, the thumb section can be easily trimmed with scissors to achieve the preferred length.

- **Breathable and low profile**
- **Adjustable support**
- **Comfort for extended daily use**

Indications: Osteoarthritis, CMC- or MCP-joint instability, inflammations, sprains, rheumatoid arthritis, after thumb arthroplasty or after surgery for thumb arthrosis, or hyper mobility in the thumb or wrist.

■ CG0108 – Neptune™ C-GO Thumb – Right / Left					Wrist circumference
Size	1 · XS	2 · S/M	3 · L/XL	4 · XXL	
cm	13–15	15–17	17–20	20–23	

Design protected in the EU under 15085791



NRX® C-GO Kids Thumb

CGO208

Providing effective thumb support for children can be challenging. During the development of the C-GO® Kids orthosis, our focus was to create an off-the-shelf solution specifically designed for children aged 4 to 12, with proportions carefully adapted to fit both pediatric thumb and wrist anatomy.



CGO208 NRX® C-GO Kids Thumb

Based on the CGO208 orthosis design, the NRX® C-GO Kids is manufactured from perforated NRX® material with a soft plush inner lining for enhanced comfort during daily wear. The outer surface features loop textile, allowing hook fasteners to attach securely anywhere on the orthosis for a fully adjustable and customized fit.

The adjustable thumb sling strap is designed to individually support and offload the CMC joint. Positioned over the thenar eminence, the sling features an inner hook surface that attaches directly to the orthosis material, enabling customized support, stabilization, and positioning of the CMC joint according to the child's needs.

NRX C-GO® Kids is available in two sizes and in left or right versions.

Indications: Thumb and thenar muscle hypoplasia, Brachial plexus injury, Cerebral palsy, Juvenile idiopathic arthritis, sprains, or hyper mobility in the thumb or wrist.



NRX®



CGO208 – NRX® C-GO Kids Thumb – Right / Left

Wrist circumference

Size	00 · XXXS kids	0 · XXS kids	
cm	9–11	11–13	

Design protected in the EU under 15085791

Orthoses are our passion

mediroyal

Mediroyal Nordic AB

Box 7052 · Staffans väg 6B · SE-192 07 Sollentuna, Sweden
Tel 08-506 766 00 · info@mediroyal.se · www.mediroyal.se

Mediroyal Asia Pacific Co. Ltd.

Unit A10, 20F., No. 282, Shizheng N. 2nd Rd., Xitun Dist., Taichung City 407, Taiwan (R.O.C.)
+886-4-2254-8069 · sales.mediroyal@gmail.com · www.mediroyal.se/en